SALTER[®]

SINCE 1760

Rechargeable Analyser Scale

Instruction Manual

EN | Please retain instructions for future reference.

Safety Instructions

This product is not a toy.

This product contains no user serviceable parts. Only a qualified electrician should attempt repairs.

Do not place, store or use the product near direct sunlight or direct heat sources or in a damp/ high humidity environment.

This product is a high precision device; handle it with care.

This product is intended for domestic use only.

Battery Safety

The supply terminals are not to be short-circuited.

In the event of a battery and cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.

Secondary cells and batteries need to be charged before use. After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.

Keep batteries out of the reach of children.

Do not crush, puncture, dismantle or otherwise damage batteries.

Do not store above 60 °C.

Do not dispose of exhausted batteries in household waste.

Do not expose the batteries to heat or fire as they may explode.

Once the batteries have been exhausted, dispose of the product at a local recycling station. Contact the local government offices for further details.

This product is fitted with a 3.7 V 1200 mAh Lithium-ion rechargeable battery (non-removable).

Care and Maintenance

Wipe the exterior of the scale with a soft, damp cloth and allow it to dry thoroughly.

Never immerse the scale in water or any other liquid. Do not use harsh or abrasive cleaning detergents or scourers to clean the scale, as this could cause damage.

Before First Use

Switch on the scale using the power switch on the bottom edge of the scale.

Position the scale on a firm flat surface.

For use on carpet attach enclosed carpet feet.

Setting the Body Analyser Scale to Zero

Position the scale on a flat, level surface. Press the centre of the platform and make sure no weight is on the scale. [0.00/ZERO] will be displayed. To select a weight unit (kg, st or lb) press [➤] or [◀]. The scale will switch off and is now ready for use.

NOTE: THE ZERO SETTING PROCESS MUST BE REPEATED IF THE SCALE IS MOVED. AT ALL OTHER TIMES STEP STRAIGHT ONTO THE SCALE.

Using the Body Analyser for a Weight Readings Only

STEP 1: Step on to the platform and stand still.

STEP 2: After 2–3 seconds, the weight reading will be displayed.

STEP 3: Step off the platform; the weight reading will be displayed for several seconds before the scale turns off.

Setting User Memory

STEP 1: Press on the platform centre with your foot and then remove.

STEP 2: Press [⊙], a user number will flash.

STEP 3: Once the user number is flashing, select by pressing [➤] or [◀]. Press [⊙] to confirm the selection.

STEP 4: The male [♂] or female [♀] symbol will flash. Press [➤] to select female [♀], male [♂], female athlete [♀♂] or male athlete [♂♂], then press [⊙].

STEP 5: The height display will flash next. Press [➤] or [◀] as necessary to set user height, then press [⊙].

STEP 6: To set user age, wait until the age display flashes. Press [➤] and [◀] to set age then press [⊙] to confirm.

STEP 7: '0.0' will be displayed. To add another user or update the memorised data, repeat steps 1–6.

ATHLETE MODE: An athlete is defined as a person who is involved in intense physical activity of approximately 12 hours per week and who has a resting heart rate of approximately 60 beats per minute or less.

NOTE: ATHLETE MODE WILL ONLY FUNCTION IF USER AGE IS SET BETWEEN 15 AND 50 YEARS.

Using the Body Analyser Scale

This scale uses BIA (Bio Impedance Analysis) technology which passes a small electrical impulse through the body to distinguish fat from lean tissue. The electrical impulse cannot be felt and is perfectly safe. Contact with the body is made via stainless steel pads on the glass platform.

NOTE: DURING PREGNANCY A WOMAN'S BODY COMPOSITION CHANGES CONSIDERABLY IN ORDER TO SUPPORT THE DEVELOPING CHILD. UNDER THESE CIRCUMSTANCES, BODY FAT PERCENTAGE READINGS COULD BE INACCURATE AND MISLEADING. PREGNANT WOMEN SHOULD THEREFORE ONLY USE THE WEIGHT FUNCTION.

For more information, visit [salter.com/understanding-your-measurements].

STEP 1: Position the scale on a flat, level surface. Press the centre of the platform and make sure no weight is on the scale.

STEP 2: Press [⊙]. Select user number by pressing [➤] or [◀].

STEP 3: The LCD display will show the saved data for that user, followed by a zero reading.

STEP 4: With bare feet, step onto the platform and stand still. Feet must be in contact with the metal pads on the platform and evenly spaced.

STEP 5: User weight, body mass index, body fat percentage [FAT %] and body water percentage [WATER %] will be displayed after a few seconds, followed by the basal metabolic rate [BMR], bone mass [BONE] and muscle mass [MUSCLE %].

NOTE: FOR CONSISTENCY, ALWAYS WEIGH ON THE SAME SCALE PLACED ON THE SAME FLOOR SURFACE.

DO NOT COMPARE WEIGHT READINGS FROM ONE SCALE TO ANOTHER AS SOME DIFFERENCES WILL EXIST DUE TO MANUFACTURING TOLERANCES. PLACING THE SCALE ON A HARD, EVEN FLOOR WILL ENSURE THE GREATEST ACCURACY AND CONSISTENCY.

IT IS RECOMMENDED THAT READINGS ARE TAKEN AT THE SAME TIME OF DAY, PREFERABLY EARLY EVENING AND BEFORE A MEAL FOR THE MOST CONSISTENT RESULTS.

THE SCALE ROUNDS UP OR DOWN TO THE NEAREST INCREMENT. IF WEIGHT IS MEASURED TWICE AND THE READINGS DIFFER BY ONE INCREMENT, THE WEIGHT LIES BETWEEN THE TWO READINGS.

Charging the Battery

When the battery is low, 'Lo' will appear on the display.

Connect the Type-C connector on the USB charging cable into the 5 V  charging port on the scale and the other end to a compatible USB port.

Whilst charging, '- - - -' will flash on the screen. '- - - -' will stop flashing once the battery has reached full charge.

Troubleshooting

'Err2' is showing on the LCD display when taking body fat and body water measurements.

The user does not have bare feet. Remove shoes and socks before using the scale. The condition of the user's skin is affecting the reading. Wipe feet with a damp cloth, leaving them slightly moist before stepping onto the scale.

Warning indicators

Lo – Charge battery.

0-Ld – Weight exceeds maximum capacity.

Err – Unstable weight. Stand still.

Err2 – Out of range, incorrect operation or poor feet contact.

Specifications

Product code: SA00507

Display: LCD

Battery: 1 x rechargeable 3.7 V Li-ion battery

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Enheter: kg/st/lb

EN If this product does not reach you in an acceptable condition, please contact our Customer Services Department at **salter.com**. Please have your delivery note to hand as details from it will be required. If you wish to return this product, please return it to the retailer from where it was purchased with your receipt (subject to their terms and conditions). To be eligible for the extended guarantee, go to **guarantee.upgs.com/salter** and register your product within 30 days of purchase.

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©Salter trademark. All rights reserved.
BRAND MERCHANT PTY LTD. SUITE 8, 8A ST ANDREWS STREET,